



GLP / GMP Conform documentation of your results.

It provides easily readable records on normal paper which last for years.

It's universal interface allows you to attach "CITIZEN" Balances & Scales via RS232C Interface.

Complete documentation of all data including calibration, weighing results, correct date and time information even after power failure.

Weighing results can be statistically evaluated. Mean value and standard deviation as well as the total weight are automatically printed out due to built-in statistics and totalization function.

CPR02
Date, time, statistic and totalization functions
Dot Matrix Printer
24 Character
Faster than 1 line per second
Normal Paper Roll, Wood Free Paper,
Width: 57.5 ± 0.5mm, Core Diameter : 8mm Min.
AC Input 100V - 300V & 50Hz
DC Output 13V & 1.5 Amps
RS 232 C
Speed of transmission: 1200, 2400, 4800, 9600, 19200 baud.
Data Bit - 8, Stop Bit 1, Parity - None, Space, Odd, Even
Working temperature : 5°C to 50°C, relative humidity: 10% - 80% (no condensation)
88.5 (H) x 135 (W) x 210 (L)
9 pin Rs232 Cable, 13V DC Adapter, Ink Ribbon
Paper Roll, user manual and warranty

WORLDWIDEMARKETING: CITIZEN SCALE (I) PVT. LTD.

Unit No. E2, Plot No. 11, WICEL Estate, Opp. Seepz Gate No.1, Andheri (E), Mumbai 400 093. (INDIA) Tel.: +91-22-4243 7700 • Fax: +91-22-4243 7800 • E:mail: sales@citizenscales.com • Web: www.citizenscales.com